

Ambu® Disposable Face Mask

DATASHEET

Description

The Disposable Face Mask is an anatomically correct shaped face mask designed for oxygenating and ventilating any patient.

The face mask comes in a full range of 6 different sizes ensuring that the user can always find a mask that fits any patient.

Features

It is a single patient use face mask that comes both with and without a check valve. The face mask has a thumb rest for easy handling while providing a tight seal.

The mask shape ensures that the inner shape of the dome is spacious enough for the mouth and nose also when pressure is applied to the mask.



Ambu Disposable Facemasks connected to Ambu SPUR II.



KEY BENEFITS

Crystal clear dome for easy observation of the patient's condition

Dome with thumb rest for an easy grip

Soft, shaped, airfilled cuff allows a tight fit to the patient's face

Available with or without hook ring; hook ring can easily be removed if not needed

The single patient use concept offers minimized risk of cross infection

All sizes are delivered individually packed, non-sterile, in packs of 20 units

Phthalate-free material

Ambu

ORDERING SPECIFICATIONS WITHOUT HOOK RING

Item no.	Type
000252051	Ambu Disposable Face Mask Neonate
000252052	Ambu Disposable Face Mask Infant
000252053H	Ambu Disposable Face Mask Toddler
000252054H	Ambu Disposable Face Mask Small adult/Child
000252055H	Ambu Disposable Face Mask Medium adult
000252056H	Ambu Disposable Face Mask Large adult
000252951	Ambu Disposable Face Mask w. check valve Neonate
000252952	Ambu Disposable Face Mask w. check valve Infant
000252953H	Ambu Disposable Face Mask w. check valve Toddler
000252954H	Ambu Disposable Face Mask w. check valve Small adult/Child
000252955H	Ambu Disposable Face Mask w. check valve Medium adult
000252956H	Ambu Disposable Face Mask w. check valve Large adult

WITH HOOK RING

Item no.	Type
000252053	Ambu Disposable Face Mask Toddler
000252054	Ambu Disposable Face Mask Small adult/Child
000252055	Ambu Disposable Face Mask Medium adult
000252056	Ambu Disposable Face Mask Large adult
000252953	Ambu Disposable Face Mask w. check valve Toddler
000252954	Ambu Disposable Face Mask w. check valve Small adult/Child
000252955	Ambu Disposable Face Mask w. check valve Medium adult
000252956	Ambu Disposable Face Mask w. check valve Large adult

SIZE

Type	Weight
Neonate	9 g
Infant	14 g
Toddler	25 g
Small adult/Child	33 g
Medium adult	47 g
Large adult	49 g

MATERIALS

Product	
Dome	Polyvinylchloride
Cuff	Polyvinylchloride
Check valve	Polycarbonate
Packaging	
Pouch	PP (Polypropylene)
Box	Cardboard
Environment	
Phthalate-free	Yes
Not made with natural rubber latex	Yes



Ambu Disposable Facemasks with check valve.



Ambu Disposable Facemasks without check valve.

Ambu

Ambu A/S

Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
ambu.com