

DATASHEET

The Sani-Child Manikin is economical and allows the training of the BLS algorithm for children.

It features realistic anatomical landmarks including sternum, rib cage and substernal notch.

It allows an optimal training of chest compression in combination with realistic ventilation. The airways open with a correct head tilt

Sanitary System

For preventing cross-contamination, a face shield/lung system is inserted easy and quickly into the manikin. The system enables a realistic exhalation and cleaning or disinfecting is no longer required.

Sani-Child Manikin



KEY BENEFITS

- Economical child manikin for BLS training
- Anatomical landmarks
- Realistic training of chest compression and ventilation
- Visible chest rise and expiration
- Airways open only in correct head position
- Sanitary system prevents cross-contamination: single-use face shield/lung system

ORDERING SPECIFICATIONS

Item no.

100-2140U	Sani-Child Manikin Torso
100-2141U	Face Shield/Lung System for Sani-Child (set of 100)

Ambu

🏢 Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com